


Please read these instructions/warnings carefully before use and keep for future reference.


**Description of this Thermometer**

1. LCD
2. On/Off button
3. Measuring tip

**Intended Use**

This digital thermometer is to be used only for measuring temperatures orally, rectally or from the armpit. It is suitable for home use.

The thermometer should be operated by persons 11+ years.

**How To Use The Thermometer**

Prior to use, check the device, ensuring that the probe tip is not damaged.

**Turning On The Thermometer**

To turn on the thermometer, press the button until the LCD screen display turns on.

**Using The Thermometer**

Once the °C symbol appears on the LCD display (after self check), place the thermometer using one of the chosen measuring methods noted below. The thermometer will read for about 30 seconds and beep four times once complete.

**Measuring Methods**
**1. Oral Measurement**

Turn on the thermometer, then place the flexible tip under the tongue (firmly against tissue). Close the mouth and breathe through the nose while the measurement is being taken. The thermometer will "Beep" four times to indicate the measurement has been taken. Remove the thermometer and read the temperature displayed on the display. **Note:** the oral temperature is 0.3°C – 0.8°C lower than the rectal temperature.


**2. Armpit Measurement**

Turn on the thermometer, place the flexible tip in the armpit parallel to the arm and press the arm tight against the body. Please ensure direct contact with skin. Hold the thermometer in this position until you hear the "Beep". The measurement has been taken and the temperature will be displayed on the LCD display. **Note:** the armpit temperature can vary up 1.5°C from than the rectal temperature.


**3. Rectal Measurement**

This is the most reliable method of measurement and is especially suitable for babies and small children. Turn on the thermometer, carefully insert the flexible tip 2cm-2.5cm into the rectum. To take the rectal temperature of a new-born baby, lay the baby on its back on a hard surface and lift its legs to an angle of approximately 90°. Hold them by the ankles in order to prevent their body from turning while taking the temperature. Remember that the internal body temperature is always slightly higher than the external temperature (about 0.5°C). You can extend the reading time by up to 1 minute if required. Wait for the thermometer to "Beep" to indicate that the measurement has been taken, then remove the thermometer and read the temperature on the LCD display temperature.


**Warnings & Limitations**

- Before first use and after every successive use, clean and sanitise the thermometer using 70% Isopropyl alcohol.
- Ensure the battery cover is in place prior to operation.
- Allow the thermometer to adapt to a new testing environment prior to use.
- This device is not intended for use on premature babies or users suffering from ulcerations, significant scarring or open wounds present at the measurement sites.
- Rectal use not recommended for those suffering from haemorrhoids.
- The thermometer is not intended for use in an oxygen rich environment.
- This device is not a toy; store in a safe place and keep out of the reach of children.
- Never leave children alone while taking their temperature.
- The thermometer can rapidly detect body temperature ranging from 32°C-42.9°C (tolerance 0.1°C - 0.2°C). If you are worried about your temperature reading please consult your Doctor.
- Do not expose the thermometer to naked flames or hot water.
- Do not drop the thermometer.
- Do not use immediately after consuming hot or cold liquids – this can cause an inaccurate temperature reading.
- Do not use the thermometer immediately after taking medication – the body can react to the medication causing the temperature to change, therefore, giving erroneous results.
- Do not use immediately after exercising, swimming or bathing.
- When measuring temperature orally, keep the mouth closed if possible to get an accurate reading. Do not remove the thermometer until the "Beep" has sounded.
- You must wait until the "Beep" has sounded before reading the measurement to ensure an accurate reading.
- Do not modify this device
- Strong electromagnetic fields can potentially affect the correct functioning of the thermometer. Do not measure while talking on the phone or while near high frequency equipment (i.e. microwave ovens etc.).

**Cleaning & Disinfecting**

The thermometer can be cleaned by wiping it with 70% Isopropyl alcohol.

In order to keep it clean and hygienic, please disinfect the thermometer before and after each use.

Do not use boiling water or abrasive cleaners to disinfect the thermometer, as this can cause damage (**Note:** The device is not intended to be sterilised).

**Calibration & Maintenance**

The digital thermometer has been calibrated during the manufacturing process and periodic calibration is not required. However, the thermometer may require testing after device impact (e.g. dropped) or there is an indication that the device accuracy limits are not as specified. If the device is damaged please dispose according to the "disposal" information found on this page.

Do not attempt any repairs to the device. You may contact the manufacturer for calibration/repairs.

**Battery Replacement**

When the battery symbol on LCD displays, please replace the battery by: Remove the battery cover. Pull out the inner frame carefully. Remove the battery with tweezers.

Only use a replacement battery of the same size and voltage (LR41 1.5V). Replace the cover.

To prolong the battery life, turn off the thermometer as soon as the reading has been taken. The thermometer will turn off automatically after approximately 3 minutes. If the device is not going to be used for a long period, we suggest removing the battery.

**Storage**

The thermometer must be stored at room temperature in a dry place, away from heat sources and direct sunlight.

If the thermometer is not being used for a prolonged period of time, please remove the battery to extend the battery life.

**Disposal**


Dispose of the device in accordance with local regulations. Dispose of at designated public collection points in EU Countries - 2012/19/EC WEEE Directive. If you have any queries, please refer to your local waste disposal authority. To protect the environment, dispose of used batteries at your local retail store or at designated public collection sites.

**Symbols**

	Follow Operating Instructions		Lot Number
	BF Type Applied Part		Date of Manufacture
	Disposal in accordance with Directive 2012/19/		Manufacturer
	IP code of the device; the grade of waterproof		Catalogue #

**Technical Data**

Storage/transport temperature: (-25°C to 55°C)

In Use temperature: (5°C to 35°C)

Relative Humidity: (<95 RH)

Measurement Range: (32.0°C to 42.9°C)

Accuracy: 0.1°C (35.5°C to 42.0°C)

0.2°C (<35.5°C, >42.0°C)

Scale: (0.1°C)

Battery Type: (LR41 – 1.5V Button cell battery)

IP22: Protection classification for the ingress of water (EN60529)

**Guarantee**

Your thermometer has a one year warranty from the date you purchase the product.

This does not include the battery. Failures resulting from the following misuse are not covered by the warranty:

- Unauthorised opening of the product.
- Damage to the product.
- Improper use of the thermometer not in accordance with these directions.

**Manufacture**

**EXCLUSIVELY  
PRODUCED FOR  
& DISTRIBUTED BY**

**CIGA Healthcare Ltd**  
Kilcran House,  
Kildowney Road,  
Ballymena, BT44 9EY, UK

**UK Customer Care Line:**  
0800 0430318  
**R.O.I. Customer Care Line:**  
0818 333 181

**Email**  
info@suresign.com  
**Online**  
www.suresign.com

**Hangzhou Hus'an Medical &  
Health Instruments Co.,Ltd.**  
Building 2,1#Fuzhu Nan RD,  
Wuchang Town, Yuhang District,  
Hangzhou, Zhejiang 310023 China

**Shanghai International Holding  
Corp. GmbH (Europe)**  
Eiffelstrasse 20, 20537 Hamburg,  
Germany